# GC Fuji PLUS CAPSULE

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### RECOMMENDED INDICATIONS




1. Cementation of metal-based inlays, onlays, crowns and bridges
2. Cementation of resin inlays, onlays, crowns and bridges
3. Cementation of all ceramic inlays
4. Cementation of high strength (zirconia based) all ceramic crowns and bridges
5. Cementation of metal, ceramic and fiber posts

### CONTRAINDICATIONS

1. Pulp capping.
2. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

### DIRECTIONS FOR USE

|                                                                 |             |
|-----------------------------------------------------------------|-------------|
| Powder / Liquid ratio (g/g)                                     | 0,36 / 0,18 |
| Mixing time (sec.)                                              | 10"         |
| Working time (min., sec.) (at 23 °C) (from start of mixing)     | 2'00"       |
| Final finishing commencing time (after seating the restoration) | 4'15"       |

1. TOOTH PREPARATION
  - a) Prepare tooth in usual manner. For pulp capping, use calcium hydroxide.
  - b) Clean the prepared teeth with pumice and water. Note : Apply GC Fuji PLUS CONDITIONER for 20 seconds using a cotton pellet or sponge.
  - c) Rinse thoroughly with water. Remove excess moisture by blotting with a cotton pellet or gently bowing with an oil-free air. DO NOT DESICCATE. Best results are obtained when prepared surfaces appear moist (glistening).
2. RESTORATION PREPARATION
 

Make sure that the restoration is pretreated and handled according to the manufacturer's instructions.
3. CAPSULE ACTIVATION AND MIXING
  - a) Before activation, shake the capsule or tap its side on a hard surface to loosen the powder (Fig. 1).
  - b) To activate the capsule, push the plunger until it is flush with the main body (Fig. 2).
  - c) Immediately place the capsule into a GC CAPSULE APPLIER and click the lever once (Fig. 3). The capsule is now activated. Note : The capsule should be activated just before mixing and used immediately.
  - d) Immediately remove the capsule and set it into a mixer (or an amalgamator) and mix for 10 seconds (+/- 4.000RPM) (Fig. 4). Note: When using a rotational mixing device, it is recommended not to use the centrifuge option.
4. DISPENSING
  - a) Immediately remove the mixed capsule from the mixer and load it into the GC CAPSULE APPLIER.
  - b) Make two clicks to prime the capsule (Fig. 5). To adjust the direction of the nozzle, hold the applicator with the capsule towards you and turn the capsule body.
  - c) Extrude directly into the restoration (Fig. 6).
  - d) To remove the used capsule, push the applicator release button. Twist the capsule and pull upwards.
5. CEMENTATION
  - a) Coat the internal surface of the restoration with sufficient cement and seat immediately. The working time is 2 minutes from start of mixing at 23°C (73.4°F). Higher temperatures will shorten working time. Note: For more complex procedures requiring longer working time, it is recommended to refrigerate the capsules in the original aluminium foil.
  - b) Maintain moderate pressure.
  - c) Start removing excess cement when the excess cement feels rubbery (Fig. 7).
  - d) Finishing can be started 4 minutes 15 seconds after seating the restoration.

**STORAGE**  
Store the capsules in the original aluminium foil in a cool and dark place 4-25°C (39,2-77,0°F). (Shelf life : 2 years from date of manufacture)

**SHADES**  
Yellow  
A3 (Vita® shade)  
Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

**PACKAGE**  
50 capsules.  
Average contents per capsule : 0.36g powder and 0.18g (0.16mL) liquid. Minimum net volume of mixed cement per capsule : 0.19mL.  
Option : GC CAPSULE APPLIER (1 piece).

- CAUTION**
1. In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water. To avoid contact, a rubber dam and/or cocoa butter can be used to isolate the operation field from the skin or oral tissue.
  2. In case of contact with eyes, flush immediately with water and seek medical attention.
  3. GC Fuji PLUS CONDITIONER contains 2% Ferric Chloride. It is not recommended to use GC Fuji PLUS CONDITIONER in combination with products containing Epinephrine as this may stain oral tissue.
  4. Do not use GC Fuji PLUS CAPSULE in combination with eugenol containing materials as eugenol may hinder GC Fuji PLUS CAPSULE from setting or bonding properly.
  5. This product is not indicated for filling or core build-up.

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Prior to use, carefully read the instructions for use.

EN

## GC Fuji PLUS CAPSULE

RADIOPAQUE REINFORCED GLASS IONOMER LUTING CEMENT IN CAPSULES

For use only by a dental professional in the recommended indications.

### RECOMMENDED INDICATIONS

- Cementation of metal-based inlays, onlays, crowns and bridges
- Cementation of resin inlays, onlays, crowns and bridges
- Cementation of all ceramic inlays, onlays, coronis and mostwose.
- Cementation of high strength (zirconia based) all ceramic crowns and bridges
- Cementation of metal, ceramic and fiber posts

### CONTRAINDICATIONS

- Pulp capping.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

### DIRECTIONS FOR USE

|                                                                 |             |
|-----------------------------------------------------------------|-------------|
| Powder / Liquid ratio (g/g)                                     | 0,36 / 0,18 |
| Mixing time (sec.)                                              | 10"         |
| Working time (min., sec.) (at 23°C) (from start of mixing)      | 2'00"       |
| Final finishing commencing time (after seating the restoration) | 4'15"       |

Fig. 1

### 1. TOOTH PREPARATION

- Prepare tooth in usual manner. For pulp capping, use calcium hydroxide liner.
- Clean the prepared teeth with pumice and water.
- Note: GC Fuji PLUS CONDITIONER for 20 seconds using a cotton pellet or sponge. Rinse thoroughly with water. Remove excess moisture by blotting with a cotton pellet or ceramic tile. DO NOT RUB. DO NOT DESICCATE. Best results are obtained when prepared surfaces appear moist (glistening).


Fig. 3

- RESTORATION PREPARATION Make sure that the restoration is pretreated and handled according to the manufacturer's instructions.

Fig. 4

- CAPSULE ACTIVATION AND MIXING a) Before activation, shake the capsule or tap its side on a hard surface to loosen the powder (Fig. 1).

Fig. 5

- Activate the capsule, push the plunger until it fits flush with the main body (Fig. 2).
- Immediately place the capsule into a GC CAPSULE APPLIER and click the lever once (Fig. 3). The capsule is now activated. Note: The capsule should be activated just before mixed and used immediately.

Fig. 6

- Immediately remove the capsule and set it into a mixer (or an amalgamator) and mix for 10 seconds (+/- 4,000RPM) (Fig. 4).

Fig. 7

- When using a rotational mixing device, it is recommended not to use the centrifuge option.

Fig. 8

- DISPENSING a) Immediately remove the mixed capsule from the mixer and load it into the GC CAPSULE APPLIER.


Fig. 9

- Make two clicks to prime the capsule (Fig. 5).
- Adjust the direction of the nozzle, hold the applicator with the capsule towards you and turn the capsule body.
- Extrude directly into the restoration (Fig. 6).
- To remove the used capsule, push the applicator release button. Twist the capsule and pull upwards.


Fig. 10

- CEMENTATION a) Coat the internal surface of the restoration with sufficient cement and seat immediately. The working time is 2 minutes from start of mixing at 23°C (73.4°F). Higher temperatures will shorten working time.


Fig. 11

- For more complex procedures requiring longer working time, it is recommended to refrigerate the capsules in the original aluminium foil.
- Maintain moderate pressure.
- Start removing excess cement when the excess cement feels rubbery (Fig. 7).


Fig. 12

- Finishing can be started 4 minutes 15 seconds after seating the restoration.


Fig. 13

- STORAGE Store the capsules in the original aluminium foil in a cool and dark place 4-25°C (39.2-77.0°F). (Shelf life: 2 years from date of manufacture)


Fig. 14

- SHADES Yellow A3 (Vita® shade) Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

Fig. 15

- PACKAGE 50 capsules Average contents per capsule : 0.36g powder and 0.18g (0.16mL) liquid. Minimum net volume of mixed cement per capsule : 0.19mL. Option : GC CAPSULE APPLIER (1 piece).

Fig. 16

- CAUTION 1. In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water. To avoid contact, a rubber dam and/or color barrier can be used to isolate the operation field from the skin or oral tissue.


Fig. 17

- In case of contact with eyes, flush immediately with water and seek medical attention.
- GC Fuji PLUS CONDITIONER contains 2% Ferric Chloride. It is not recommended to use GC Fuji PLUS CONDITIONER in combination with products containing Epinephrine as this may stain oral tissue.
- Do not use GC Fuji PLUS CAPSULE in combination with eugenol containing materials as eugenol may hinder GC Fuji PLUS CAPSULE from setting or bonding properly.
- This product is not indicated for filling or core build-up.


Fig. 18

- CEMENTATION a) Coat the internal surface of the restoration with sufficient cement and seat immediately. The working time is 2 minutes from start of mixing at 23°C (73.4°F). Higher temperatures will shorten working time.


Fig. 19

- For more complex procedures requiring longer working time, it is recommended to refrigerate the capsules in the original aluminium foil.
- Maintain moderate pressure.
- Start removing excess cement when the excess cement feels rubbery (Fig. 7).


Fig. 20

- STORAGE Store the capsules in the original aluminium foil in a cool and dark place 4-25°C (39.2-77.0°F). (Shelf life: 2 years from date of manufacture)


Fig. 21

- SHADES Yellow A3 (Vita® shade) Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.


Fig. 22

- PACKAGE 50 capsules Average contents per capsule : 0.36g powder and 0.18g (0.16mL) liquid. Minimum net volume of mixed cement per capsule : 0.19mL. Option : GC CAPSULE APPLIER (1 piece).


Fig. 23

- CAUTION 1. In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water. To avoid contact, a rubber dam and/or color barrier can be used to isolate the operation field from the skin or oral tissue.

Fig. 24

- In case of contact with eyes, flush immediately with water and seek medical attention.
- GC Fuji PLUS CONDITIONER contains 2% Ferric Chloride. It is not recommended to use GC Fuji PLUS CONDITIONER in combination with products containing Epinephrine as this may stain oral tissue.
- Do not use GC Fuji PLUS CAPSULE in combination with eugenol containing materials as eugenol may hinder GC Fuji PLUS CAPSULE from setting or bonding properly.
- This product is not indicated for filling or core build-up.

Fig. 25

- CEMENTATION a) Coat the internal surface of the restoration with sufficient cement and seat immediately. The working time is 2 minutes from start of mixing at 23°C (73.4°F). Higher temperatures will shorten working time.


Fig. 26

- For more complex procedures requiring longer working time, it is recommended to refrigerate the capsules in the original aluminium foil.
- Maintain moderate pressure.
- Start removing excess cement when the excess cement feels rubbery (Fig. 7).

Fig. 27

- STORAGE Store the capsules in the original aluminium foil in a cool and dark place 4-25°C (39.2-77.0°F). (Shelf life: 2 years from date of manufacture)

Fig. 28

- SHADES Yellow A3 (Vita® shade) Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

Fig. 29

- PACKAGE 50 capsules Average contents per capsule : 0.36g powder and 0.18g (0.16mL) liquid. Minimum net volume of mixed cement per capsule : 0.19mL. Option : GC CAPSULE APPLIER (1 piece).

Fig. 30

- CAUTION 1. In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water. To avoid contact, a rubber dam and/or color barrier can be used to isolate the operation field from the skin or oral tissue.

Fig. 31

- In case of contact with eyes, flush immediately with water and seek medical attention.
- GC Fuji PLUS CONDITIONER contains 2% Ferric Chloride. It is not recommended to use GC Fuji PLUS CONDITIONER in combination with products containing Epinephrine as this may stain oral tissue.
- Do not use GC Fuji PLUS CAPSULE in combination with eugenol containing materials as eugenol may hinder GC Fuji PLUS CAPSULE from setting or bonding properly.
- This product is not indicated for filling or core build-up.

Fig. 32

- CEMENTATION a) Coat the internal surface of the restoration with sufficient cement and seat immediately. The working time is 2 minutes from start of mixing at 23°C (73.4°F). Higher temperatures will shorten working time.

Fig. 33

- For more complex procedures requiring longer working time, it is recommended to refrigerate the capsules in the original aluminium foil.
- Maintain moderate pressure.
- Start removing excess cement when the excess cement feels rubbery (Fig. 7).

Fig. 34

- STORAGE Store the capsules in the original aluminium foil in a cool and dark place 4-25°C (39.2-77.0°F). (Shelf life: 2 years from date of manufacture)

Fig. 35

- SHADES Yellow A3 (Vita® shade) Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

Fig. 36

- PACKAGE 50 capsules Average contents per capsule : 0.36g powder and 0.18g (0.16mL) liquid. Minimum net volume of mixed cement per capsule : 0.19mL. Option : GC CAPSULE APPLIER (1 piece).

Fig. 37

- CAUTION 1. In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water. To avoid contact, a rubber dam and/or color barrier can be used to isolate the operation field from the skin or oral tissue.

Fig. 38

- In case of contact with eyes, flush immediately with water and seek medical attention.
- GC Fuji PLUS CONDITIONER contains 2% Ferric Chloride. It is not recommended to use GC Fuji PLUS CONDITIONER in combination with products containing Epinephrine as this may stain oral tissue.
- Do not use GC Fuji PLUS CAPSULE in combination with eugenol containing materials as eugenol may hinder GC Fuji PLUS CAPSULE from setting or bonding properly.
- This product is not indicated for filling or core build-up.

Fig. 39

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