Prior to use, carefully read the instructions for use.

EN

GC Fuji II LC CAPSULE

RADIOPAQUE LIGHT CURED REINFORCED GLASS IONOMER RESTORATIVE IN CAPSULES

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

- Restoration of Class III, V and limited Class I cavities.
- Restoration of primary teeth.
- Core build up.
- Cases where a radiopaque restoration is required.
- Geriatric applications.
- As a base or liner.

CONTRAINDICATIONS

- Pulp capping.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

DIRECTIONS FOR USE

Powder / Liquid Ratio (g/g)	0,33 / 0,10
Mixing Time (sec.)	10"
Working Time (min., sec.)	3'15"
Light Curing Time (sec.)	20"
Depth of Cure (A2) (mm)	1,8

Test conditions : Temperature (23 +/-1°C)
Relative humidity (50 +/-10%)
ISO 9917-2 : 1998 (E) (Light-activated cements) (Type I)

1. CAVITY PREPARATION

- Prepare tooth using standard techniques. Extensive mechanical retention is not necessary. For pulp capping, use calcium hydroxide.
- Wash and dry but DO NOT DESICCATE. GC CAVITY CONDITIONER or GC DENTIN CONDITIONER is recommended to remove the smear layer (Fig. 1).
- CAPSULE ACTIVATION AND MIXING**
 - Before activation, shake the capsule or tap its side on a hard surface to loosen the powder (Fig. 2).
 - To activate the capsule, push the plunger until it flush with main body (Fig. 3).
 - Immediately place the capsule into a metal GC Capsule Applier and click the lever once (Fig. 4). The capsule is now activated. Note : The capsule should be activated just before mixing and used immediately.
 - Immediately remove the capsule and set it into a mixer (an amalgamator) and mix for 10 seconds at high speed (+/-4.000 RPM) (Fig. 5).

3. RESTORATIVE TECHNIQUE

- Immediately remove the mixed capsule from the mixer and load it into the GC Capsule Applier.
- Make two clicks to prime the capsule then syringe (Fig. 6).
The working time is 3 minute 15 seconds from start of mixing at 23°C (73.4°F). Higher temperature will shorten working time. Note :
 - To adjust the direction of the nozzle, hold the applier with the capsule towards you and turn the capsule body.
 - To remove the used capsule, push the applier release button. Twist the capsule and pull upwards.
- Remove surface moisture but DO NOT DESICCATE.
- Extrude cement directly into preparation (Fig.7). Avoid air bubbles.
- Form the contour and place a matrix if required.
- Light-cure for 20 seconds using a suitable visible light curing device (470nm wavelength) (Fig. 8). Place light source as closely as possible to the cement surface. Note : For cavities deeper than 1.8 mm, use a layering technique.
- FINISHING**
Finish under water spray using superfine diamond bur, silicone point and polishing strips (Fig. 9). Note : Apply a final coat of GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure for 10 sec.) or G-COAT PLUS (light cure for 20 sec.)(Fig. 10).

STORAGE

Store in a cool and dark place (4-25°C / 39.2-77.0°F). (Shelf life : 2 years from date of manufacture)

SHADE

A1, A2, A3, A3.5, A4, B2, B3, B4, C2, C4, D2 (11 shades based on Vita® shades).

Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

PACKAGES

- Single shade package : 50 capsules.
- Assortment package : 50 capsules (10 pcs each of A4, B2, B3, C2; 5 pcs each of B4, D2).
*Average contents per capsule : 0.33g powder and 0.10g (0.085mL) liquid.
*Minimum net volume of mixed cement per capsule : 0.10mL.
Option : GC Capsule Applier (1 piece).

CAUTION

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water.
- In case of contact with eyes, flush immediately with water and seek medical attention.

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CE 0086

GC

MANUFACTURED BY
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

DISTRIBUTED BY
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

EU: GC EUROPE N.V.
Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33,
B-3001 Leuven, Belgium TEL: +32 16 74 10 00

GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
TEL: +1-708-597-0900

GC ASIA DENTAL PTE. LTD.
19 Loyang Way, #06-27 Singapore 508724
TEL: +65 6545 7588

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