Prior to use, carefully read the instructions for use.

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GC Fuji II

RADIOPAQUE GLASS IONOMER RESTORATIVE CEMENT

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

- Restoration of primary teeth.
- Core build-up.
- Restorations of Class III, V and limited Class I cavities.

CONTRAINDICATIONS

- Pulp capping.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

DIRECTIONS FOR USE

- POWDER AND LIQUID DISPENSING**
 - The standard powder to liquid ratio is 2.7/1.0g, 1 level scoop of powder to 1 drop of liquid.
 - For accurate dispensing of powder, lightly tap the bottle against the hand. Do not shake or invert.
 - Hold the liquid bottle vertically and squeeze gently.
 - Cure bottles immediately after use.
- MIXING**
 - Dispense powder and liquid onto the pad. Using the plastic spatula, divide the powder into 2 equal parts. Mix the first portion with all of the liquid for 10 seconds. Incorporate the remaining portion and mix the whole thoroughly for 15-20 seconds (total time within 30 seconds).
- RESTORATIVE TECHNIQUE**
 - Prepare tooth using standard techniques. Extensive mechanical retention is not necessary. For pulp capping use calcium hydroxide.
 - Wash and dry but DO NOT DESICCATE.
 - GC CAVITY CONDITIONER** or **GC DENTIN CONDITIONER** is recommended to remove the smear layer.
 - Mix the required amount of cement. Working time is 1 minute 45 secs at 23°C (73.4°F). Higher temperatures will shorten working time.
 - Remove surface moisture, but DO NOT DESICCATE.
 - Transfer cement to the preparation using a syringe or suitable placement instrument. Avoid air bubbles.
 - Form the contour, and place a matrix if required during the first 2 minutes of setting.
 - Setting time is 5*30" from start of mixing. Once set, remove the matrix and immediately apply GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure).
- FINISHING**
 - One visit treatment
 - After initial set, correct the contour WITHOUT WATER SPRAY using standard techniques.
 - Apply GC Fuji VARNISH or GC Fuji COAT LC and wait for about 15 minutes.
 - Final finishing and polishing under water spray can now be carried out.
 - Apply a final coat of GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure for 10 sec.) or G-COAT PLUS (light cure for 20 sec.).
 - Standard method
 - 24 hours after placement finish and polish the restoration under water spray using standard techniques.

STORAGE
Store at temperature of 4-25°C (39.2-77.0°F). (Shelf life : 3 years from date of manufacture).

SHADE

No. 21 - pale yellow
No. 22 - yellow brown
No. 23 - dark grey
No. 26 - deep brown
M - Modifier

PACKAGES

1. Bottle of 15g powder with measure.
2. Bottle of 10g liquid (8.0 mL).
3. 1-1 pkg. : 15g powder, 10g (8.0mL) liquid, powder measure, mixing pad (No. 22).
ISO 9917 (Dental water based cements)

CAUTION

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- DO NOT mix powder or liquid with any other glass ionomer components.

Last revised : 01/2009

MANUFACTURED BY
GC CORPORATION
76-1 Hasunuma-cho, Habashi-ku, Tokyo 174-8585, Japan

DISTRIBUTED BY
GC CORPORATION
76-1 Hasunuma-cho, Habashi-ku, Tokyo 174-8585, Japan

EU: GC EUROPE N.V.
Researchpark Haasrode-Leuven 1240, Interfeulvenlaan 33,
B-3001 Leuven, Belgium TEL: +32 16 74 10 00

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