




G-aenial™ Bond

UNIVERSAL BONDING AGENT

1 TRULY UNIVERSAL ONE-STEP BONDING AGENT



Now available
in unit dose!



TOUCH of GENIUS™

- Top rated 7th generation, one-step self-etch bonding agent
- Improved bond strengths over other leading bonding agents to both enamel and dentin
- Choose from selective-etch, self-etch or total-etch
- Easy application - 30 seconds from start to finish

Dental PRODUCT SHOPPER
RECOMMENDED
2011

THE DENTAL ADVISOR
2012
Top 7th-generation Bonding Agent

THE DENTAL ADVISOR
Editors' Choice

THE DENTAL ADVISOR
2012 PREFERRED PRODUCTS

Higher Bond Strength, Ultimate Versatility

Economical 300 Dropper Bottle

G-ænial™ Bond's small nozzle design minimizes evaporation of solvent. Cost per drop is \$0.41. Compare to other competitors:


G-ænial™ Bond	\$0.41
Xeno IV†	\$0.65
Prime & Bond NT†	\$0.70
Clearfil S3 Bond†	\$0.70
Brush & Bond†	\$0.70
iBond Total Etch†	\$0.92

Based on 2014 retail pricing of bottle refills per mL.




See and Feel the Difference

After preparation of G-ænial™ Bond, you can see a frosty, matte-like surface effect. The first layer of composite sticks nicely to G-ænial™ Bond.




Traditional adhesives leave a wet, glossy surface after light activation. Composites do not adhere well to the adhesive, making it hard to manipulate.



G-ænial™ Bond creates a sticky, frosty, matte-like surface to aid in composite adhesion and manipulation.

G-ænial™ Bond is HEMA-Free!

The HEMA-free formulation avoids degradation of the collagen fibrils. Water is not attracted to the tooth surface, improving bond strength. The dentin/adhesive interface, the nano-interaction zone, and the presence of hydroxyapatite crystals within this zone, together with the HEMA-free formulation, are key factors for the superior long-lasting bond strength found with G-ænial™ Bond.



Use with G-ænial™ Universal Flo, the First Truly Universal Injectable Composite



004220 G-ænial™ Bond Refill
Contains: One G-ænial Bond Bottle (5 mL).

008483 G-ænial™ Bond Unit Dose
Contains: Box of 50 Unit Doses (0.1 mL) and 50 disposable applicators (fine).

004217 G-ænial™ Bond Starter Kit
Contains: One G-ænial Bond Bottle (5 mL), 50 disposable applicators (fine) and 20 disposable dishes.

¹Data on file.
†Not a registered trademark of GC America Inc.

