

ADDENT, INC.

The AdDent, Inc. Polarizing Filter is attached to the Rite-Lite 2 High CRI with four magnets.

The Polarizing Filter can be purchased with the Rite-Lite 2 High CRI or factory installed on some older models. Check with factory.

Parts List:

P/N 110075 - Rite-Lite 2 High CRI
w/o Polarizing Filter

P/N 110075-PF - Rite-Lite 2 High CRI
w/Polarizing Filter

P/N 640083 - Polarizing Filter

P/N 640015 - Neutral Color Pad (Pk. of 25)

What Doctors and Ceramists say:

*"Great product. I have used it to verify my own shade selections. Every single shade selection I chose was verified correct by Rite*lite 2. This would be an indispensable tool for anyone who doesn't have a great eye for color."*

---Stephen Balloch, DDS
Glastonbury, CT

"We love this light. We haven't been off on one shade since we started using this light."

---Dennis Lempke
Lempke Dental Laboratories
Maplewood, MN

*"Using the Rite*lite 2 made me change my color choices at least 30% of the time."*

---David Zalucky, DDS
Wilton, CT

**Rite-Lite 2 Tri-Spectra/Hi CRI
Shade Matching Light**

Features:	Benefits:
3 light settings for shade taking and color analysis	Use daylight setting to acquire shade; use room and ambient light settings to verify shade
Hi Color Rendition Index (CRI)	Recreates visual spectrum like that of natural sunlight
Polarizing Filter	Enhances visualization of internal tooth structure
3 lighting conditions from one source	Easily toggle between 3 modes with one-touch control
Standardizes lighting conditions between multiple locations	Provides consistent lighting conditions between dental office and dental lab for shade analysis
Controlled viewing area	Eliminates problems of shadowing from upper lip and distracting outside color sources from interfering with shade selection
Neutralizing color tabs included	Eliminates color distraction in the oral cavity
Battery operated	No warm up time required
L.E.D. provides full power instantly	L.E.D. provides 30,000 hours of life
Voltage Regulator	Constant light output with low battery indicator

Rite*lite2™
TRI-SPECTRA/Hi CRI SHADE LIGHT



*Clearly, a Bright Idea.
Now with Hi CRI!*



43 MIRY BROOK RD. DANBURY, CT 06810
203-778-0200 • WWW.ADDENT.COM

4/2015 PN640132

