



FREEGENOL™ Temporary Pack - Temporary Cement

PRODUCT DESCRIPTION

FREEGENOL is a patented non-eugenol temporary cement with ideal characteristics for all types of provisional restorations. This temporary cement is completely compatible and safe to use with all acrylic resin materials as well as with plastic temporary restorations. It is bio-compatible with tooth structure and soft tissue and will not deteriorate or compromise the integrity of temporary restorations. FREEGENOL temporary cement is easy to use and mixes to a smooth, creamy consistency. Due to a highly thixotropic characteristic, the cement stays where it is placed, yet flows well under pressure. FREEGENOL sets hard and fast in only three minutes and allows removal of provisionals to be accomplished quickly and effortlessly.



FEATURES AND BENEFITS

Patented eugenol-free formula

- Will not burn or irritate soft tissue.
- Acrylic restorations will not soften or deteriorate.
- No cracking or crazing of plastic temporaries.

Thixotropic properties

- Cement stays where placed; flows well under pressure.

Easy mix paste

- Mixes to a smooth, creamy consistency in 20-30 seconds.

Sets hard and fast in 3 minutes

- Saves time.
- Allows fast, easy removal of flashing.
- Set time can be adjusted by modifying amount of accelerator.

Easy removal of temporaries

- Minimal time and effort required to prepare for final restoration.

PRODUCT CODE

136501 FREEGENOL Temporary Cement - Standard Package

Contains: One tube Base (50gm), one tube Accelerator (20gm), mixing pad



GC AMERICA INC.

3737 WEST 127TH STREET
ALSIP, IL 60803
TEL. (800) 323-7063
FAX (708) 371-5103